



## **✓ Applications**

The range of compressors is designed to test, adjust and calibrate measuring instruments mechanical or electronic pressure by comparison (pressure gauges, pressure transmitters or sensors or switches).

These generators are equipped with a pressure tank visible level, two valves connected by a set of two-stage variable volumes. The capstan can generate the pressure pushing the fluid through the variable volumes to fittings provided for receiving the pressure instruments.

These compressors are designed for industry or laboratory. They are robust and easy to use and have a very high long-term stability.

They have the opportunity to grow in dead weight tester by simply adding a kit with : a set piezometric and a set of weights.

### ✓ Technical specifications of the range CH1

- Fluid : Mineral oil or synthetic oil depending on the model colorless compatible with

medical or food uses

**- Tank** : 125 cm3

- Volume de fluid expelled : 48 cm<sup>3</sup>

- Check/adjust the seat : leveling bubble and adjustable feet

- Connection for the instrument to be tested : swivel G1/2 standard – other optional fittings

- Carter: Light alloy aluminium AG3 + high-resistance paint

Capstan: Removable for transport
 Operating temperature: 18 to 28°C – 64 to 82°F

#### ✓ Technical product specifications

- **Model** : CH1-1250B

- **Measuring range**: 0 to 1250 bar / 0 to 20000 PSI / 0 to 125000 kPa

- **Sensitivity**: 10% Pmaxi = 0.1 bar / 50% Pmaxi = 0.2 bar / 100% Pmaxi = 0.5 bar

- Weight comparator : 18 k

- Fluid: Oil H40 - fluid density: 860 kg/m<sup>3</sup>


- **Section of piston**: low pressure = 1.1 cm<sup>2</sup>

- Number of piston :

Material Piston/Cylinder:
 P = stainless steel / C = treated steel
 Piping:
 copper+ polyamide + treated steel

# ✓ Divers

- Manuel: a detailed operating instructions supplied with the device + Certificate + 0.5 liter Oil
- Maintenance: Our technical department is at your disposition for the revision, calibration or service your unit Calibration device recommended every 2 years.



## **✓ Options of the generator**

- Kit for dead weight (CH2 -> BH2)
- Pump hands with tank (Ref. OP0158)
- Pressure microajusteur : increases the accuracy (contact us)
- adaptation corrosives fluid (skydroll, brake fluid)
- Oil H15: 1 liter (ref.OP0013) or 5 liters (ref.OP0010)
- Suitcase for the geenrator (ref.OP0002):









505 à 615 505 to 615

✓ Accessories (ask us our specific documentation)

- Connectors up to 1200 bar: suitcase connectors M (ref.OP0174) - suitcase connectors G (ref.OP0171) - suitcase connectors NPT (ref.OP0172) - suitcase connectors BSP-TR (ref.OP0173) - suitcase with 17 connectors M + G + NPT + BSP-TR (ref.OP0037) - unit connector







# ✓ Transport and packaging

- Packaging : pack is provide for the shipping

Designation/Reference	Dimension / carton or shipping crate	Weight empty / total weight (packaging + device)	
Comparator	470 x 470 x 240 mm	Carton empty 1.5 kg Total weight ±20 kg	
packaging SB0001	480 x 480 x 340 mm	Box empty 10 kg Total weight ±30 kg (comparator + packaging)	



SB0001

- Note : shipping and taxes are your responsibility - EX-Works Vendôme (France)

## √ Other models available in the range CH1

(Datasheet is available on request or on our website : www.aremeca-instrumentation.com)

		Measuring range		
Models	Accuracy	Bar or kg/cm²	PSI	kPa
CH1-125B	10 <sup>-3</sup>	0 to 125	0 to 1800	0 to 12500